*Samples mounted in
ClaroFast*

*Proben eingebettet in
ClaroFast*

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Table of Contents	Page	English
Hot Mounting Materials.....	5	
Mount Release Agent	5	
Accessories	5	
Optimizing the Process Time.....	6	
Helpful Hints	7	
Mounting Data.....	8	
Trouble-Shooting.....	10	

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Hot Mounting Materials

Product	Application	Specific property	Material	Contents	Cat. no.
ClaroFast	Glass clear mounts Porous samples Surface electrical insulator for ConduFast	Transparent Medium shrinkage	Acrylic Thermoplastic	1 kg 7.5 kg 25 kg	40100055 40100054 40100053
CitoFast	Fast mounting times For soft materials	Very fast mounting. Low process times also when it is used as 'backing' for DuroFast or LevoFast	Acrylic resin with aluminium filler Thermoplastic Grey	1 kg 7.5 kg	40100068 40100069
ConduFast	Electrolytic polishing	Electrically conductive Very low shrinkage	Acrylic resin with iron filler Thermoplastic Black/Grey	1 kg	40100039
DuroFast	Edge retention Planeness	Very low shrinkage Good adhesion Very low removal rate	Epoxy with mineral and glass filler Thermosetting Black	1 kg 7.5 kg	40100044 40100045
LevoFast	Edge-retention for soft to medium hard materials	Very low shrinkage High removal rate	Melamine with mineral and glass filler Thermosetting Light yellow	1 kg 7.5 kg	40100057 40100058
PolyFast	Edge retention	Very low shrinkage Medium removal rate	Phenolic resin with carbon filler Thermosetting Black	1 kg 7.5 kg	40100036 40100037
	Examination in scanning electron microscope	Electrically conductive Low emission in the vacuum chamber			
PuriFast	Universal use Backing material	Medium shrinkage High removal rate	Polypropylene material with wood flour and mineral filler Thermoplastic Uncolored	2.5 kg 7.5 kg 30 kg 90 kg	40100082 40100083 40100084 40100085
PuriFast Green	Universal use Backing material Color marking	Medium shrinkage High removal rate	Polypropylene material with wood flour and mineral filler Thermoplastic Green	2.5 kg 7.5 kg 30 kg 90 kg	40100086 40100087 40100088 40100089
PuriFast Red	Universal use Backing material Color marking	Medium shrinkage High removal rate	Polypropylene material with wood flour and mineral filler Thermoplastic Red	2.5 kg 7.5 kg 30 kg 90 kg	40100090 40100091 40100092 40100093
MultiFast Black	Routine examination Backing material	Medium shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting Black	2.5 kg 7.5 kg 25 kg 75 kg	40100064 40100065 40100066 40100067
MultiFast Green	Routine examination Color marking	Medium shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting Green	2.5 kg 7.5 kg 25 kg 75 kg	40100078 40100079 40100080 40100081
MultiFast Red	Routine examination Color marking	Medium shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting Red	2.5 kg 7.5 kg 25 kg 75 kg	40100074 40100075 40100076 40100077

Mount Release Agent

Product	Application	Material	Contents	Cat. no.
AntiStick	To prevent the mounts sticking to the rams	Very fine stearate powder	8 g	40300043

Note. Do not use AntiStick with PuriFast

Accessories


Product	Application	Material	Contents	Cat. no.
Fixation Clips 9 mm	Clip for holding small and thin specimens when mounting	Metal spring	100 pcs.	40300026
Fixation Clips 6 mm		Metal spring	100 pcs.	40300025
EpoDye		Fluorescent dye	20 g	40300002
EasyDoser	Manual dosing unit for hot mounting materials	For Struers 1 kg containers		05796902

Optimizing the Process Time


When specimens are to be mounted in series it is recommended to optimize the heating and cooling times.

One of the factors limiting both the heating and cooling times is the relatively low heat conductivity of the resins. An efficient way to reduce the process time is to minimize the distance the heat needs to travel through the resin. When mounting metallic pieces (with high heat conductivity), the heating and cooling times may be reduced when:


- The amount of resin is optimized:


- A relatively high instead of a low sample is mounted:

- Choose smallest possible mounting cylinder:

- Good thermal contact between metal and ram is ensured:

When optimizing the time it should be reduced stepwise with intermediate inspection of the mount. Insufficient heating and cooling times will result in artifacts (See: "Trouble Shooting").

Helpful Hints

ConduFast and ClaroFast

For electrolytical purposes, ClaroFast, although not electrically conductive, can be used in conjunction with ConduFast. By adding a small amount of ClaroFast first, and then completing the required amount with ConduFast, a mount with a conducting body and an insulating preparation surface is formed.

DuroFast and LevoFast with PuriFast, MultiFast and CitoFast

DuroFast or LevoFast can be used with PuriFast, MultiFast or CitoFast as backing material. Add approx. 1/3 of DuroFast/LevoFast to make the preparation surface and then complete the mount with 2/3 of the backing material. Use CitoFast for fast mounting.

Mount Release Agent

This is recommended to be applied to the mounting rams as a thin layer before the mounting process begins. This prevents the adhesive qualities of the resins from making it difficult to remove the mounts afterwards. Mount release agent, AntiStick cannot be used with PuriFast.

Distance to Cylinder Wall

The distance between the specimen and the cylinder wall must be a minimum of 3 mm (1/8"), to avoid cracks in the resin. This is especially critical for specimens with sharp corners.

Small Specimens

Small, thin specimens can be supported during the hot mounting process by the use of Struers Fixation clips. Only metal clips should be used for hot mounting.

Clean Specimens

To obtain the best results the specimens must be clean, dry and free from grease. If necessary, clean with alcohol or another suitable degreasing fluid.

Preheating

For porous and/or pressure sensitive specimens, such as minerals, electronic parts etc., it is advantageous to soften the mounting material by heating, before applying pressure. Preheating is also useful when using thermoplastic materials i.e. ClaroFast and PuriFast.

We recommend applying pressure for the last minute only to protect the specimen. If the total heating time is i.e. 8 minutes, preheating should be set to 7 minutes at 0 bar and hereafter 1 minute at 100 - 350 bar according to the properties of the specimens.

Temperature Sensitive Specimens

The temperature for the mounting process can be reduced to a minimum of 150°C for all mounting materials. This is useful when dealing with temperature sensitive materials. For PuriFast the heating time is unchanged using 150 °C or 180 °C.

For very temperature sensitive specimens, hot compression mounting should be avoided. Use Struers cold mounting resins instead.

Porous Specimens

Thermoplastic mounting material, (ClaroFast and PuriFast), penetrate into porous specimens. The best results are obtained by initially preheating the material (see Preheating).

For very porous materials, hot compression mounting should be avoided. Use of Struers epoxy cold mounting resins is recommended.

Pressure Sensitive Specimens


Thermoplastic materials, (ClaroFast and PuriFast), should be used. The best results are obtained by initially preheating the materials (see Pre-heating).

For very pressure sensitive specimens, hot compression mounting should be avoided. Use Struers cold mounting materials instead.

Specimen Removal

It is crucial that the mounting material and parameters are correctly matched to each individual specimen. If needed ClaroFast and PuriFast can be reprocessed.

Specimens mounted in thermosetting mounting material cannot be reprocessed. They have to be re-mounted after removing the mounting material using a hammer or bench vice, which might damage the specimen.

*Polished zinc coating,
Mounted in PolyFast
1000 x*

Mounting Data

The recommended mounting parameters appear in the following pages.

CitoPress-5/-15/-30

The heating and cooling times in the tables refer to the following conditions:

- The heating times in the tables refer to the full process time counted from the start of the process and not from when the pre-set temperature is reached. Add ½ - 1 minute to the heating time for the first mount to heat the cylinder.
- The mounting parameters are based on specimens with a volume approx. 20% of the total volume of the mount. If larger specimens are mounted, less resin should be used.
- If smaller specimens or specimens with low heat conductivity are mounted, the heating and cooling times should be increased. It might also be necessary to increase the pressure to avoid pores in the cured mount.
- A higher pressure is recommended with ClaroFast to avoid "cotton ball" effect.

- The amount of mounting material in the tables has been adjusted to result in a final height of the mounts of approx. 20 mm (0.8").
- If several mounting materials are combined in one mount, use the process parameters for the mounting materials with the longest times. When using CitoFast as 'backing', use the CitoFast process parameters. For several, or very complicated specimens, add one min. of heating time.

On all CitoPress models, three preset cooling levels are available: High, Medium and Low.

Cylinder dia.	Resin		Heating			Cooling		Time
25 mm / 1"	Type	Quantity	Time	Temp.	Pressure	Time	Rate	Total time
		[ml]	[min]	[°C]	[bar]	[min]		[min]
	ClaroFast	20	4	180	350	6	Low	9.5
	CitoFast	15	2	180	300	1	High	3
	ConduFast	15	3	180	300	1	High	4
	DuroFast	15	3	180	325	2	High	5
	LevoFast	20	3	180	350	1	High	4
	PolyFast	15	3.5	180	325	1.5	High	5
	PuriFast	20	3	150 - 180	300	2	High	5
	MultiFast	20	3	180	300	2	High	5

Cylinder dia.	Resin		Heating			Cooling		Time
30 mm	Type	Quantity	Time	Temp.	Pressure	Time	Rate	Total time
		[ml]	[min]	[°C]	[bar]	[min]		[min]
	ClaroFast	20	4	180	350	6.5	Low	10.5
	CitoFast	25	2.5	180	300	1	High	3.5
	ConduFast	20	3.5	180	250	1.5	High	5
	DuroFast	20	3.5	180	325	2	High	5.5
	LevoFast	25	3.5	180	250	2	High	5.5
	PolyFast	20	3.5	180	250	1.5	High	5
	PuriFast	25	3	150 - 180	250	2	High	5
	MultiFast	25	3	180	250	2	High	5









Cylinder dia.	Resin		Heating			Cooling		Time
1¼"	Type	Quantity	Time	Temp.	Pressure	Time	Rate	Total time
		[ml]	[min]	[°C]	[bar]	[min]		[min]
	ClaroFast	20	4	180	350	6.5	Low	10.5
	CitoFast	25	2.5	180	300	1	High	3.5
	ConduFast	20	3.5	180	250	1.5	High	5
	DuroFast	20	4	180	325	2	High	6
	LevoFast	25	3.5	180	250	2	High	5.5
	PolyFast	20	3.5	180	250	1.5	High	5
	PuriFast	25	3.5	150 - 180	250	2	High	5.5
	MultiFast	25	3.5	180	250	2	High	5.5

Cylinder dia.	Resin		Heating			Cooling		Time
1½"	Type	Quantity	Time	Temp.	Pressure	Time	Rate	Total time
		[ml]	[min]	[°C]	[bar]	[min]		[min]
	ClaroFast	30	4	180	350	7	Low	11
	CitoFast	45	3	180	300	1.5	High	4.5
	ConduFast	35	3.5	180	250	2	High	5.5
	DuroFast	35	4.5	180	350	2.5	High	7
	LevoFast	40	4.5	180	250	2	High	6.5
	PolyFast	30	4.5	180	250	2	High	6.5
	PuriFast	40	4.5	150 - 180	250	2.5	High	7
	MultiFast	40	4.5	180	250	2.5	High	7







Cylinder dia.	Resin		Heating			Cooling		Time
40 mm	Type	Quantity	Time	Temp.	Pressure	Time	Rate	Total time
		[ml]	[min]	[°C]	[bar]	[min]		[min]
	ClaroFast	35	4	180	350	6.5	Low	10.5
	CitoFast	45	3	180	300	1.5	High	4.5
	ConduFast	40	3.5	180	250	2	High	5.5
	DuroFast	40	4.5	180	350	2.5	High	7
	LevoFast	50	5	180	250	2	High	7
	PolyFast	35	4	180	250	2	High	6
	PuriFast	45	4	150 - 180	250	2.5	High	6.5
	MultiFast	45	4	180	250	2.5	High	6.5


Cylinder dia.	Resin		Heating			Cooling		Time
50 mm / 2"	Type	Quantity	Time	Temp.	Pressure	Time	Rate	Total time
		[ml]	[min]	[°C]	[bar]	[min]		[min]
	ClaroFast	55	5	180	250	8.5	Low	13.5
	CitoFast	65	3.5	180	250	2	High	5.5
	ConduFast	65	4	180	250	2	High	6
	DuroFast	70	5.5	180	250	3	High	8.5
	LevoFast	75	6	180	250	3	High	9
	PolyFast	55	5.5	180	250	2	High	7.5
	PuriFast	70	6	150 - 180	250	4	High	10
	MultiFast	70	6	180	250	4	High	10

Trouble-Shooting

General Problems					
Radial Cracking 	Cause: Insufficient distance between specimen edge/corner and cylinder wall, or specimen has sharp corners Solution: Increase cylinder diameter or reduce specimen size. The distance between the specimen and the cylinder wall must be a minimum of 3 mm to avoid cracks in the resin. This is especially critical for specimens with sharp corners.	Shrinkage 	Cause: Incorrect choice of resin. Solution: Re-mount a new specimen using a resin with a lower linear shrinkage value.	Blistering 	Cause: Insufficient heating time. Solution: Increase heating time, or increase process temperature. Cause: Overcured surface. Solution: Decrease process temperature. Cause: Entrapped gas within mount. Solution: Preheat resin.
Bulging 	Cause: Insufficient cooling. Solution: Increase cooling time.	Porosity 	Cause: Excessive temperature. Solution: Reduce process temperature.	Voiding within large mounts 	Cause: Insufficient heating time. Solution: Increase heating time. Cause: Excessive temperature. Solution: Reduce process temperature. Cause: Insufficient force/pressure. Increase mounting force/pressure.
Dull Surface Finish 	Cause: Insufficient heating time. Solution: Increase heating time.	Adhesion between mount and rams 	Cause: Insufficient application of mould release agent. Except for PuriFast. Release agent cannot be used with PuriFast. Solution: Apply mould release agent. This must always be applied to the mounting rams as a thin layer before the mounting process begins. This prevents the resins from sticking to the rams and makes it easier to remove the mounts afterwards. Cause: Insufficient heating time. Solution: Increase heating time. Cause: Excessive force/pressure. Solution: Decrease mounting force/pressure.	Individual Grains Visible on Mount* 	Cause: Resin has cured without force/pressure. Solution: Increase force/pressure during heating cycle. Cause: Insufficient heating time. Solution: Increase the heating time and/or temperature. *Thermosetting resins only ** For PuriFast are the individual grains always visible.

Trouble-Shooting

ConduFast Problems			
Lack of conductivity 	Cause: No contact with specimen due to use of an excessive amount of ClaroFast, (refer to Helpful Hints). Solution: Re-mount a new specimen using a smaller amount of ClaroFast than before. Cause: Insufficient heating time Solution: Reinsert the mount in the press and reprocess using an increased heating time.	Metallic particles in resin removed during the electrolytic process 	Cause: Insufficient amount of ClaroFast, (refer to Helpful Hints). Solution: Re-mount a new specimen using a larger amount of ClaroFast. Cause: Excessive grinding time. Solution: Re-mount with a new specimen.
ClaroFast Problems			
Internal cracking 	Cause: Excessive cooling rate. Solution: Reduce the cooling rate	“Cottonball” effect in centre of mount 	Cause: Excessive cooling rate. Solution: Reduce the cooling rate. Cause: Insufficient heating time. Solution: Reduce the physical height of the mount, lower the process temperature and increase the heating time. Cause: Moist resin. Solution: Dry the resin by exposing the open container to 30-70°C for 2 hours.
LevoFast Problems			
The mounts turn dark after preparation 	Cause: Insufficient heating time Solution: Increase heating time and/or temperature		
PolyFast and MultiFast Problems			
Parts of mount become light or colourless on contact with alcohol 	Cause: Insufficient heating time. Solution: Increase heating time and/or temperature		
PuriFast problem			
Cracking specimens when clamping in specimen holder.	Cause: Specimen clamped with too much force. Solution: Use max. 1.5 N or tighten the screw until it touches the mount and then turn the Allen key 2 sides (120 °) Make attention not to have a sharp edge on the specimen where the screw is clamped into the mount. If a PuriFast mount cracks it can be reprocessed to obtain a perfect mount.		

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Struers Application Guide

The following restrictions should be observed, as any violation of the restrictions may cause cancellation of Struers legal obligations:

Struers Application Guide may only be used in connection with Struers products and equipment covered by the Application Guide.

Struers assumes no responsibility for errors in the Application Guide text/illustrations. The information in the Application Guide is subject to changes without notice.

The Application Guide may mention consumables, accessories or parts not included in the present version of the Struers product you are using.

The contents of this Application Guide are the property of Struers. Reproduction of any part of the Application Guide without the written permission of Struers is not allowed.

All rights reserved. © 10th revised edition, Struers 2011.

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