

	DP-Produkt	Verfügbare Korngrößen
	DP-Suspension P	15 µm, 9 µm, 6 µm, 3 µm, 1 µm, 0,25 µm, 0,1 µm
	DP-Suspension M	15 µm, 9 µm, 6 µm, 3 µm, 1 µm
	DP-Suspension A	15 µm, 9 µm, 6 µm, 3 µm, 1 µm, 0,25 µm
	DP-Spray P	35 µm, 15 µm, 9 µm, 6 µm, 3 µm, 1 µm, 0,25 µm
	DP-Spray M	15 µm, 9 µm, 6 µm, 3 µm, 1 µm
	DP-Stick P	15 µm, 9 µm, 6 µm, 3 µm, 1 µm, 0,25 µm
	DP-Paste P	15 µm, 9 µm, 6 µm, 3 µm, 1 µm, 0,25 µm
	DP-Paste M	15 µm, 9 µm, 6 µm, 3 µm, 1 µm, 0,25 µm



TECHNISCHE DATEN

Diamantpolierprodukte (DP-Produkte)

Diamanten	Auf Wasser- oder Alkoholbasis	Gebindegröße*
Polykristallin	Wasser	125 ml** 500 ml 2,5 l 5 l
Monokristallin	Wasser	500 ml 2,5 l 5 l
Polykristallin/ Monokristallin	Alkohol	500 ml
Polykristallin	Alkohol	150 ml
Monokristallin	Alkohol	150 ml
Polykristallin	Fest	24 g
Polykristallin	Fest	10 g
Monokristallin	Fest	10 g

* Nicht alle μm -Größen sind in allen Gebindegrößen erhältlich.
Alle Einzelheiten finden Sie im Verbrauchsmaterialkatalog von Struers

** In einer Pumpflasche zur manuellen Anwendung