

Preparation disc with cone		Cat. no.
MD-Disc, 200 mm dia.		06086401
Magnetic fixation discs enable you to change the preparation surface fast and easy. It furthermore ensures less resilience when preparing your specimens. Only for use with our MD surfaces. With cone for easy replacement.		
Aluminum disc, 200 mm dia.		06086404
For use when preparing your specimens with OP and DP surfaces and other self-adhesive grinding and polishing consumables. With cone for easy replacement.		
Wet grinding disc, 200 mm dia.		06206924
The wet grinding disc is for fixation of SiC paper.		
Splash guard*		Cat. no.
For semi-automatic preparation with a MD-Disc or an Aluminium disc		06206923
For semi-automatic preparation on LaboPol-20. LaboPol comes standard with guard for manual preparation, with small gap to prevent stuck specimens		
For Wet grinding disc		06206926
For use with ø200 mm Wet grinding disc. LaboPol comes standard with guard for manual preparation, with small gap to prevent stuck specimens		
Lid for preparation disc		06206912
For protection of LaboPol-20 disc during non-use		
Disposable bowl liner 5 pcs.		49900061
For easy cleaning of your bowl		

* LaboPol-20 comes standard with a splash guard for manual preparation, with a small gap to prevent stuck specimens

Suitable for storage of 9 MD preparation surfaces in 300 mm (12"), 250 mm (10") or 200 mm (8") dia.



Storage cabinet for 7 MD preparation surfaces in 200 mm (8") dia.

