

Models	Cat. no.
Tegramin-30	06036127
Tegramin-30 with cover	06036227
Tegramin-30 with safety cover	06036527


Tegramin dosing module	Cat. no.
With 1 DP pump	06036900

Easy connector (s) included.



With 2 DP pump	06036901
----------------	----------

Easy connector(s) included.



With 2 DP pumps and 1 OP pump	06036902
-------------------------------	----------

Easy connector(s) included.



Bottle tray for 8 bottles	06036911
----------------------------------	----------

For organization of bottles.
Bottle tray with space for eight ½ or 1 liter square bottles.



Preparation discs		Cat. no.
MD-Disc, 300 mm dia.		06086403
Magnetic fixation discs enable you to change the preparation surface fast and easy. It furthermore ensures less resilience when preparing your specimens. Only for use with our MD surfaces.		
Aluminum disc, 300 mm dia.		06086406
For use when preparing your specimens with OP and DP surfaces and other self-adhesive grinding and polishing consumables.		
Cast iron disc. 300 mm dia.		06206917
Concentrically grooved, cone-type disc for lapping mineralogical specimens. Made of special SiC-resistant cast iron alloy		
Cover		06036903
Transparent cover		
Lid for preparation disc		06036906
For protection of preparation disc during non-use		
Splash guard for manual preparation		06036905
With small gap to prevent stuck specimens		

Disposable bowl liner 5 pcs.	49900052
------------------------------	----------

For easy cleaning of your bowl



Shift valve	06036904
-------------	----------

To be used if you have a recirculation system and use a method with OP-suspensions



Storage cabinet	Cat. no.
-----------------	----------

MD-Concert	05306101
------------	----------

Suitable for storage of 9 MD preparation surfaces in 300 mm (12"), 250 mm (10") or 200 mm (8") dia.



Levelling device	Cat. no.
Uniforce levelling device without pressure foot Levelling device for levelling of specimens in specimen holders 140, 160 and 200 mm dia. Without pressure foot for positioning of specimens	07046101
Uniforce levelling device with pressure foot Levelling device for levelling of specimens in specimen holders 140, 160 and 200 mm dia. With pressure foot for positioning of specimens	07046102

Specimen mover plate					Cat. no.
Coupling					06086901
For specimen mover plates.					
					
Hole size	Holes	Shape	Photo	Material	Cat. no.
25 mm	6			Stainless steel	06086932
30 mm	6			Stainless steel	06086933
40 mm	6			Stainless steel	06086934
50 mm	6			Stainless steel	06086935
1 "	6			Stainless steel	06086939
1 ¼ "	6			Stainless steel	06086940
1 ½ "	6			Stainless steel	06086941
2 "	6			Stainless steel	06086942
20x30 mm	6			Stainless steel	06086944
Accustop-30	3			Stainless steel	06086937
Accustop-40 and TYNDS holder	3			Stainless steel	06086938
-	0	-		Stainless steel	06086936

Specimen holders					
Hole size	Holes	Shape	Photo	Material	Cat. no.
25-40 mm	6			Stainless steel	06066909
25-40 mm	6			Aluminum	06066912
10-25,5 mm	12			Stainless steel	02606917
10-25 x 34 mm	6			Stainless steel	02606918
20-40 x 68 mm	3			Stainless steel	02606920
Max. 24 x 52 mm	3	-		Stainless steel	02606924
29-32mm	8			Aluminum	02606954
38-40 mm	6			Aluminum	02606955
50-51 mm	4			Stainless steel	02606963
-	0	-		Stainless steel	02606926
-	0	-		Aluminum	02606928
Flexible specimen holder					Cat. no.
For unmounted odd-shaped specimens.					07066961
Area: Min. 20 x 40 mm. Max. 100 x 100 mm Height: Min. 10 mm. Max. 20 mm Special software required. Tegramin software version ≥ 1.33					
Viaholder					Cat. no.
Preparation of PCB coupons					05726903
Adjustable precision specimen holder with diamonds stops. For up to 6 x6 PCB coupons					