

Models	Cat. no.
Tegramin-25	06026127
	
Tegramin-25 with cover	06026227
	
Tegramin-25 with safety cover	06026527
	


Tegramin dosing module	Cat. no.
With 1 DP pump	06036900

Easy connector (s) included.



With 2 DP pump	06036901
----------------	----------

Easy connector(s) included.



With 2 DP pumps and 1 OP pump	06036902
-------------------------------	----------

Easy connector(s) included.



Bottle tray for 8 bottles	06036911
----------------------------------	----------

For organization of bottles.
Bottle tray with space for eight ½ or 1 liter square bottles.



Preparation disc with cone	Cat. no.
MD-Disc, 250 mm dia.	06086402

Magnetic fixation discs enable you to change the preparation surface fast and easy. It furthermore ensures less resilience when preparing your specimens. Only for use with our MD surfaces. With cone for easy replacement.



Aluminum disc, 250 mm dia.	06086405
-----------------------------------	----------

For use when preparing your specimens with OP and DP surfaces and other self-adhesive grinding and polishing consumables. With cone for easy replacement.



Cover	06036903
--------------	----------

Transparent cover.



Lid for preparation disc	06036906
---------------------------------	----------

For protection of preparation disc during non-use.



Splash guard for preparation disc	06026905
--	----------

For manual preparation with small gap to prevent stuck specimens.



Disposable bowl liner 5 pcs.	49900052
------------------------------	----------

For easy cleaning of your bowl



Shift valve	06036904
-------------	----------

To be used if you have a recirculation system and use a method with OP-suspensions



Storage cabinet	Cat. no.
-----------------	----------

MD-Concert	05306101
------------	----------

Suitable for storage of 9 MD preparation surfaces in 300 mm (12"), 250 mm (10") or 200 mm (8") dia.



Levelling device	Cat. no.
For 140 mm dia. specimen holders	05256903

For levelling of specimens in specimen holder's 140 mm dia.



Uniforce levelling device without pressure foot	07046101
---	----------

Levelling device for levelling of specimens in specimen holders 140, 160 and 200 mm dia.
Without pressure foot for positioning of specimens



Uniforce levelling device with pressure foot	07046102
--	----------

Levelling device for levelling of specimens in specimen holders 140, 160 and 200 mm dia.
With pressure foot for positioning of specimens



Specimen mover plate					Cat. no.
Coupling					06086901
For specimen mover plates.					
					
Hole size	Holes	Shape	Photo	Material	Cat. no.
25 mm	6			Stainless steel	06086902
30 mm	6			Stainless steel	06086903
40 mm	6			Stainless steel	06086904
50 mm	3			Stainless steel	06086905
1 "	6			Stainless steel	06086909
1 ¼ "	6			Stainless steel	06086910
1 ½ "	6			Stainless steel	06086911
2 "	3			Stainless steel	06086943
Accustop 30	3			Stainless steel	06086907
-	-	-		Stainless steel	06086906

Specimen holders					
Hole size	Holes	Shape	Photo	Material	Cat. no.
10-32 mm	6	●		Stainless steel	06066917
10-32 mm	6	●		Aluminum	06066920
25-40 mm	4	●		Stainless steel	06066906
25-40 mm	4	●		Aluminum	06066903
29-32 mm	6	●		Aluminum	02606952
-	0	-		Aluminum	02606950
Flexible specimen holders					07066941
For unmounted odd-shaped specimens. Area: Min. 20 x 40 mm. Max. 90 x 90 mm Height: Min. 10 mm. Max. 20 mm Special software required. Tegramin software version ≥ 1.33					